Blom-Singer®

voice restoration systems

37-478-01 Rev. B

Indwelling Voice Prosthesis, 16 Fr.	6	ENG
Verweil-Stimmprothese, 16 Fr.	14	DEU
Prótesis de Voz Alojada, 16 Fr.	23	ESP
Prothèse Vocale à Demeure, 16 Fr.	31	FRA
Protesi Fonatoria Fissa, 16 Fr.	40	ITA
Verblijf Stemprothese, 16 Fr.	49	NED
Prótese Vocal Fixa, 16 Fr.	58	POR

INHEALTH
TECHNOLOGIES

Clinician Information Data Sheet

Informationsblatt für Ärzte/Hoja de Información
del Facultativo/Fiche technique d'information
cliniciens/Scheda informativa per il medico/
Informatieblad voor de clinicus/Folheto
Informativo do Médico

Indwelling Voice Prosthesis, 16 Fr.

Verweil-Stimmprothese / Prótesis de voz alojada / Prothèse vocale à demeure /
Protesi fonatoria fissa / Verblijf stemprothese / Prótese vocal fixa

One-way Flap Valve

Einwegeklappenventil/Válvula de aleta de una sola salida/
Clapet antiretour/Valvola a cerniera unidirezionale/
Eenrichting klepafsluiting/
Válvula de charneira unidireccional

Esophageal Retention Flange

Speiseröhrenhalterung/Brida de retención esofágica/Bride de fixation œsophagienne/Flangia di fissaggio esofageo/Slokdarm retentiefens/
Flange de retenção esofágica

Body

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




INDWELLING VOICE PROSTHESIS — CLEAN AND READY-TO-USE

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BLOM-SINGER® INDWELLING VOICE PROSTHESIS, 16 FR.

The Blom-Singer Indwelling Voice Prosthesis is intended for laryngectomees who are unable or reluctant to perform the routine removal and insertion necessary for the cleaning and maintenance of traditional patient-maintained voice prostheses. The Indwelling voice prosthesis is clinician-maintained; it is inserted and replaced on an as needed basis by a speech pathologist or physician only.

INDICATIONS FOR USE

The Blom-Singer Indwelling Voice Prosthesis is indicated for tracheoesophageal voice rehabilitation following total laryngectomy where placement or replacement of a voice prosthesis is performed by a physician or qualified clinician.

CONTRAINDICATIONS

The Indwelling Voice Prosthesis and accessories are medical products that should be used only by professionals with experience and training in their use and care. The flushing device is designed for use only in cleaning particulate matter from the voice prosthesis.

HOW SUPPLIED

The Indwelling Voice Prosthesis is available in a **Clean and ready-to-use** series for Secondary Placement. It is supplied nonsterile in a single-wrapped package and includes: Indwelling voice prosthesis, inserter, gel cap introducer, gel caps, lubricant, and a patient care kit with flushing device and instructions for use. The Flushing Device is also available separately as a 3-pack (IN4050), clean and ready-to-use.

The Blom-Singer Indwelling voice prosthesis is only sold to physicians and licensed speech pathologists and is not sold directly to the user (patient).

INSTRUCTIONS FOR USE IN SECONDARY PLACEMENT

The voice prosthesis and its accessories should always be thoroughly clean and dry to help avoid introducing contaminants into the tracheostoma, tracheoesophageal puncture, or esophagus.

The clinician's hands should always be gloved and clean and the accessories should always be clean. Removal, insertion, and cleaning of devices in the puncture should only be done with a bright light focused directly on the stoma.

Warning: Caution should be exercised to avoid accidental displacement of the device which could result in aspiration of the prosthesis into the windpipe. Should aspiration occur, the patient should attempt to cough the device out of the windpipe. Further medical attention may be necessary if coughing the device out is unsuccessful.

The following procedure is provided by Eric D. Blom, Ph.D., and Ronald C. Hamaker,

M.D., F.A.C.S., as instructions for loading the voice prosthesis into the gel cap and inserting the voice prosthesis into the tracheoesophageal puncture. Use of gloves and protective eyewear is recommended for the following procedure.

Be certain that your hands and the voice prosthesis are completely dry before beginning the following procedure. Failure to do so may cause premature dissolving of the gel cap. Always handle the voice prosthesis with clean, gloved hands. **Note:** Prior to use of the voice prosthesis, check the valve mechanism to be sure it is intact and working properly. The flap valve should close flat against the seating surface of the valve seat inside the voice prosthesis.

Please refer to the diagrams located at the front of this data sheet.

Gel Cap Loading

1. Fold the esophageal retention flange tightly in half toward the center of the valve (diagram 1).
2. Fold the two outer edges of the esophageal retention flange tightly against the folded flange (diagram 2) and hold the prosthesis in this folded position.
3. Insert the folded portion of the esophageal retention flange into the gel cap (diagram 3). Gently tuck and push the folded prosthesis as far as possible into the gel cap.
4. Use the flange introducer to gently tuck the remaining portion of the folded esophageal retention flange completely into the gel cap (diagram 4). Do not use any sharp or serrated instruments to avoid damaging the prosthesis.
5. The esophageal retention flange should be folded in a forward position and completely enclosed inside the gel cap (diagram 5).

Voice Prosthesis Insertion

1. Place the voice prosthesis on the inserter and lock the strap onto the safety peg.
2. Caution the patient not to swallow as you remove the previously inserted 18 Fr. tracheoesophageal puncture dilator to avoid aspiration of saliva into the trachea prior to inserting the voice prosthesis.
3. Apply a light coating of water-soluble lubricant (oil-free) to the tip of the gel-capped end of the voice prosthesis and immediately insert the voice prosthesis completely into the tracheoesophageal puncture by aligning the tip of the voice prosthesis partially in the puncture with the neck strap oriented upwards (diagram 6). **Note:** Do not use petroleum-based products such as Vaseline® to lubricate the prosthesis. These products can damage silicone.
4. Hold the voice prosthesis in this position of complete insertion for at least 3

minutes. This allows ample time for the gel cap to dissolve and release the esophageal retention flange within the esophagus. Patients should be instructed to swallow their saliva or warm water to facilitate dissolving the gel cap.

Caution: If the voice prosthesis does not insert easily on the first attempt, do not continue to try to insert it. Instead, reinsert a clean 18 Fr. tracheoesophageal puncture dilator for a few minutes to dilate the tracheoesophageal tract, and then try again.

5. Detach the voice prosthesis strap from the safety peg on the inserter. Place a finger against the strap and carefully withdraw the inserter from the prosthesis with a twisting motion. **Warning:** Caution should be exercised when removing the inserter to avoid accidental displacement of the voice prosthesis which could result in aspiration of the prosthesis. Secure the voice prosthesis with adhesive tape across the strap so it will not slide out.

6. Have the patient take a breath, occlude the stoma, and produce voice. Be certain there is no leakage through the prosthesis by having the patient drink water.

Optional Strap Detachment

The voice prosthesis is designed to include the optional detachment of the strap by a physician or qualified speech pathologist following confirmation that the prosthesis esophageal retention flange is securely deployed and correctly positioned within the esophageal lumen against the anterior wall of the esophagus. It is recommended to orient the voice prosthesis with the strap up.

The following procedure for confirming deployment of the esophageal retention flange is recommended by Eric D. Blom, Ph.D. and Ronald C. Hamaker, M.D., F.A.C.S.:

After insertion of the voice prosthesis in the puncture tract and completion of the 3 minute waiting period for gel cap dissolution and esophageal retention flange deployment, rotate the voice prosthesis on the inserter while it is in the puncture. If the esophageal retention flange has deployed and has properly seated against the mucosa of the anterior esophageal wall, the prosthesis will rotate easily in the puncture tract on the inserter. If the gel cap has dissolved and incorrectly deployed the esophageal retention flange inside the puncture tract rather than inside the lumen of the esophagus, considerable resistance to prosthesis rotation will be felt.

A second method for confirming gel cap dissolution and esophageal retention flange deployment is to assess the capacity to produce tracheoesophageal voice. If the voice prosthesis is correctly positioned and the esophageal retention flange deployed, airflow diverted from the trachea by digital occlusion of the tracheostoma will open the flap valve and pass up the esophagus to initiate voice.

If the proper deployment and positioning of the esophageal retention flange cannot be confirmed by one or both of the methods described above, troubleshoot gel cap dissolution and esophageal retention flange deployment by manually opening the one-way flap valve in the esophageal end of the voice prosthesis with gentle contact. Direct a bright light into the prosthesis. **Slowly** and **carefully** advance the shaft end of a long cotton swab into the voice prosthesis until **gentle contact** opens the flap valve. A flap valve that does not open with gentle contact indicates that the undeployed retention flange is preventing valve opening.

If the assessment methods described above fail to provide satisfactory verification of correct esophageal retention flange deployment, remove the voice prosthesis, dilate and remeasure the puncture tract, and repeat the voice prosthesis insertion and confirmation process.

Once deployment of the esophageal retention flange has been confirmed, detach the strap from the tracheal retention flange at the area of reduced strap width by cutting carefully with scissors (diagram 7).

Note: if you are uncertain that the esophageal retention flange has properly deployed against the anterior wall of the esophagus after using the steps described above, do not detach the strap from the prosthesis. Securely tape the strap to the neck and observe the patient for a few days for retention ability. Removal of the strap is an optional procedure. It is recommended that the strap not be detached from an indwelling voice prosthesis if the clinician anticipates that simultaneous use of a laryngectomy tube may potentially dislodge the voice prosthesis.

REMOVAL OF THE VOICE PROSTHESIS

The Blom-Singer Indwelling voice prosthesis may be removed at any time if there is any dissatisfaction with the voice prosthesis or if a different approach to speech restoration is desired. In this event, the clinician will remove the voice prosthesis and allow the tracheoesophageal puncture close permanently.

The voice prosthesis may be left in place in the tracheoesophageal puncture until it ceases to function correctly, that is, until it leaks or is not providing adequate voice for speech, or, for not longer than 6 months, whichever comes first. If the voice prosthesis is not functioning properly, the patient should return it to the clinician for evaluation.

Removal of the voice prosthesis should only be done by grasping the outer tracheal flange securely with a locking hemostat. Pull gently and firmly until the prosthesis is fully withdrawn. Insert an 18 Fr. tracheoesophageal puncture dilator and tape it in position for five minutes prior to inserting a new voice

prosthesis that has been attached to an inserter. **Never** remove one voice prosthesis and insert another voice prosthesis without first dilating the tracheoesophageal puncture with the 18 Fr. tracheoesophageal puncture dilator, and, remeasuring the tract to confirm voice prosthesis length. **Always** use a gel cap on the tip of the voice prosthesis to enable easy, atraumatic insertion. **Warning:** Clearly inform the patient **never** to attempt to remove the voice prosthesis by him/herself or allow anyone other than a qualified clinician to remove the voice prosthesis.

CLEANING AND CARE INSTRUCTIONS

To keep the inside of the voice prosthesis body and valve clean of dried phlegm, it is recommended that the prosthesis be cleaned carefully on a daily basis with the supplied flushing device. The flushing device provides a means for cleaning the voice prosthesis while in-situ, i.e., in the patient's tracheoesophageal puncture. The flushing device is designed for use only in cleaning particulate matter from the voice prosthesis.

InHealth relies on the physician or trained Speech Pathologist to instruct patients on the use of this device and to provide patients with the Blom-Singer Indwelling voice prosthesis User Information Data Sheet supplied with the product.

The following instructions should be made clear to the patient as part of the routine care and cleaning of the voice prosthesis:

Caution: Use only supplied cleaning device. Do not insert foreign objects into prosthesis. The cleaning accessories and user's hands should **always be thoroughly clean** to help avoid introducing contaminants into the tracheostoma, puncture, or esophagus. Cleaning of the voice prosthesis should only be done in front of a mirror with a bright light focused directly on the stoma so that the open end of the voice prosthesis is clearly visible. Before use, the flushing device should be rinsed with 3% hydrogen peroxide. Use long handled tweezers to carefully remove any dried phlegm that may be in the open end of the voice prosthesis.

Flushing Device

To fill the flushing device with water, squeeze the bulb and then release it in clean water. Shake excess water from the tip of the flushing device before flushing the voice prosthesis. Carefully and gently insert the tip of the flushing device into the voice prosthesis until it fits completely against the opening of the prosthesis (diagram 8). Do not flush water through the voice prosthesis until the tip of the flushing device is properly inserted into the prosthesis, thus avoiding leaking water into the trachea. Should a small amount of water leak into the trachea, this will cause coughing. Squeeze the bulb on the flushing device to flush a rapid

squirt of water through the voice prosthesis. After flushing, remove the flushing device carefully. If water will not readily squirt through the voice prosthesis, this indicates that it may be plugged with dried phlegm. Allow a few drops of water to dissolve the dried matter for a few minutes and then flush again with the flushing device until the debris is loosened.

Routine Care

Repeat cleaning regimen as needed. The flushing device should be cleaned, rinsed and air-dried after each use. Squeeze and fill the flushing device with 3% hydrogen peroxide after each use. Shake to remove any residual liquid. Allow the flushing device to dry in a well-lit, aerated area. If desired, this accessory can be washed in soap and water, rinsed and air-dried.

WARNINGS AND PRECAUTIONS

In some patients, the inability to relax the muscles of the throat may account for their inability to speak fluently and with minimal effort. This problem requires assessment by the surgeon and speech pathologist.

Patients requiring postoperative radiation may have transient interruption of voice in the third or fourth week of treatment. The voice prosthesis can be left in place or a catheter can be inserted as a stent until radiation treatments are completed.

The strap on an Indwelling voice prosthesis should only be removed after the physician or trained speech pathologist has verified the retention flange position within the esophagus. Never attempt to insert or reinsert an Indwelling voice prosthesis that has the strap removed.

Use only the specially-sized gel caps provided by InHealth Technologies. Be certain that the gel cap is completely dry and securely installed on the tip of the voice prosthesis prior to insertion into the tracheoesophageal puncture.

When using the Blom-Singer Gel Cap Insertion System, always lubricate the tip of the gel-capped end of the voice prosthesis with a water-soluble lubricant prior to insertion. Insert the voice prosthesis immediately after lubricating the gel-capped tip to prevent premature dissolving of the gel cap. Use only water-soluble lubricants (i.e. K-Y® Jelly, Surgilube® or Lubrifax®). Do not use petroleum-based lubricants, such as Vaseline®, as these can damage silicone.

Use only supplied cleaning device. Do not insert foreign objects into prosthesis.

The flushing device is designed for use only in cleaning particulate matter from the voice prosthesis. If there are tears, cracks, or structural damage to the flushing device, discontinue use.

Do not flush liquid through the voice prosthesis until the tip of the flushing device is

properly inserted into the prosthesis, thus avoiding leaking water into the trachea. Should a small amount of water leak into the trachea, this will cause coughing.

Caution should be exercised when inserting or removing cleaning devices to avoid accidental displacement of the voice prosthesis which could result in aspiration of the voice prosthesis. In the unlikely event that this should occur, the patient should attempt to cough the voice prosthesis out of the windpipe and immediately seek medical attention if unsuccessful.

If the voice prosthesis is accidentally dislodged from the puncture, the patient should be instructed to immediately place an 18 Fr. Blom-Singer tracheoesophageal puncture dilator in the puncture to keep the puncture from closing. The patient should then return to his/her clinician for reinsertion of a voice prosthesis within 24 hours.

Some patients experience the growth of a yeast (*Candida albicans*) in their mouth and/or esophagus. Yeast deposits can cause valve deformation and failure, i.e., liquid leakage through the prosthesis and an increase in the pressure necessary to voice. If a patient experiences recurrent leakage of their voice prosthesis, he/she should return to his/her clinician for evaluation. Some patients can be treated successfully for yeast with oral medication or mouth rinses.

Clearly inform the patient never to attempt to remove or reinsert the voice prosthesis under any circumstances by him/herself or allow anyone other than the clinician to do so.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician or a licensed speech pathologist.

COMPLICATIONS

Although rare, the following complications have been identified to occur with silicone voice prostheses of the Blom-Singer type. They include: stoma contamination or sepsis, which may require removal of the prosthesis and/or appropriate antibiotics; accidental aspiration of the voice prosthesis into the airway, which may require removal by a physician; aspiration of liquids through the prosthesis into the airway, which may cause the patient to cough; occasional extrusion of the voice prosthesis, requiring replacement after dilation of the tracheoesophageal puncture and additional supervision of the stoma care regimen; puncture dilation resulting in leakage of fluids around the prosthesis; inflammatory reaction around the puncture site and formation of granulation tissue; dislodgment of the voice prosthesis and subsequent closure of the tracheoesophageal puncture; intractable leakage around the voice prosthesis, requiring surgical revision or closure of the puncture; dysphagia, or difficulty swallowing; tearing or other damage to the voice prosthesis from improper use; yeast deposits causing leakage through the prosthesis or valve incompetence;

accidental ingestion of the prosthesis into the esophagus and/or GI tract.

SPECIAL ORDER PRODUCTS

If this data sheet accompanies a Special Order Product there will possibly be differences in the physical characteristics between the product enclosed and the product descriptions in this data sheet. These differences will not affect the safety or efficacy of the special order product. Special order products are nonreturnable.

ORDERING INFORMATION

USA

Blom-Singer products may be ordered directly from International Healthcare Technologies (InHealth) at the following toll free number: (800) 477-5969, Monday — Thursday, 9:30 am—7:00 pm; and Friday, 9:30 am—5:00 pm, Eastern Standard Time. The toll free number may be used from the continental USA, Alaska, Hawaii, Puerto Rico, and the Virgin Islands. Products also may be ordered at the InHealth website: www.inhealth.com

Canada

In Canada, please contact Auto Control Medical at (800) 461-0991, or Fax (800) 699-5936.

International

For a list of international distributors, please contact: International Healthcare Technologies, 1110 Mark Avenue, Carpinteria, CA 93013-2918, USA. Tel: +01(805) 684-9337, Fax: +01(805) 684-8594.

Consumer Affairs

If you have any questions or dissatisfaction with a product please contact customer service via telephone, fax, post, or e-mail. Tel: (800) 477-5969; Fax: (888) 371-1530; Post: 1110 Mark Avenue, Carpinteria, CA 93013-2918 USA; E-mail: productcomplaints@inhealth.com

RETURNED GOODS POLICY

All returned merchandise must have a Return Merchandise Authorization number (RMA), and must be unopened and undamaged. RMA numbers are obtained from InHealth Customer Service. Returns without an RMA will not be accepted. Products must be returned in unopened packages, with manufacturer's tamper evident seals intact to be accepted for replacement or credit. Products will not be accepted for replacement or credit if they have been in possession of the customer for more than 90 days. Special order products are nonreturnable. There is a 20% restocking fee on all returned merchandise.

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BIBLIOGRAPHY

1. Blom, E.D.: "Some Comments on the Escalation of Tracheoesophageal Voice Prosthesis Dimensions", **Arch Otolaryngol**, 129: 500-502, 2003.
2. Eerenstein, S.E., Grolman, W., Schouwenburg, P.F.: "Downsizing of Voice Prosthesis Diameter in Patients with Laryngectomy", **Arch Otolaryngol**, 128: 838-841, 2003.
3. Blom, E.D., Singer M.I., and Hamaker, R.C., **Tracheoesophageal Voice Restoration Following Total Laryngectomy**. San Diego: Singular Publishing Group, Inc.1998.
4. Blom, E.D., Hamaker, R.C.: "Tracheoesophageal Voice Restoration Following Total Laryngectomy," **Cancer of the Head and Neck**, E.N. and Suen, J. (Eds). W.B. Saunders Publishers. 839-852, 1996.
5. Blom, E.D., "Tracheoesophageal Speech," **Seminars in Speech and Language**. New York: Thieme Medical Publishers. 16:3, 191-204, 1995.
6. Pauloski, B.R., Blom, E.D., Logemann, J.A., Hamaker, R.C. "Functional Outcome After Surgery for Prevention of Pharyngospasms in Tracheoesophageal Speakers. Part II: Swallow Characteristics", **Laryngoscope**, 105:10: 1104-1110, 1995.
7. Blom, E.D., Pauloski, B.R., Hamaker, R.C. "Functional Outcome After Surgery for Prevention of Pharyngospasms in Tracheoesophageal Speakers. Part I: Speech Characteristics", **Laryngoscope**, 105:10: 1093-1103, 1995.
8. Blom, E.D., Hamaker, R.C., Freeman, S.B.: "Postlaryngectomy Voice Restoration." Lucente, F.E. (Ed.): **Highlights of the Instructional Courses**, Vol. 7. St. Louis: Mosby-Year Book. pp 3-10, 1994.
9. Blom, E.D., Singer M.I., and Hamaker, R.C. "A Prospective Study of Tracheoesophageal Speech", **Arch Otolaryngol**, 112: 440-447, 1986.
10. Blom, E.D., Singer, M.I. "Disinfection of Silicone Voice Prostheses", **Arch Otolaryngol**, 112: 1986.
11. Blom, E.D., Singer M.I., and Hamaker, R.C. "An Improved Esophageal Insufflation Test", **Arch Otolaryngol**, 111: 211-212, 1985.
12. Hamaker, R.C., Singer M.I., Blom E.D. and Daniels, H.A. "Primary Voice Restoration at Laryngectomy", **Arch Otolaryngol**, 111: 182-186, 1985.
13. Singer, M.I. and Blom, E.D. "A Selective Myotomy for Voice Restoration After Total Laryngectomy", **Arch Otolaryngol**, 107: 670-673, 1981.
14. Singer, M.I. and Blom, E.D. "An Endoscopic Technique for Restoration of Voice After Laryngectomy", **Ann Otol Rhinol Laryngol**, 89:6: 529-533, 1980.

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Manufacturer

Hersteller/Fabricante/Fabricant/
Produttore/Fabrikant

International Healthcare Technologies

A Division of Helix Medical Inc.
1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
www.inhealth.com

European Representative

Europäische Vertretung/Representante para
Europa/Représentant européen/Rappresentante
europeo/Europese vertegenwoordiger/
Representante Europeu

EMDAR BV

IJsselburcht 3, POB 5486
6802 EL Arnhem, The Netherlands

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Blom-Singer is a registered trademark of Hansa Medical Products, Inc., InHealth is a registered trademark of International Healthcare Technologies, A Division of Helix Medical, Inc., © 2004 International Healthcare Technologies / Blom-Singer ist ein eingetragenes Warenzeichen der Hansa Medical Products Inc., InHealth ist ein eingetragenes Warenzeichen von International Healthcare Technologies, ein Unternehmen der Helix Medical, Inc., / Blom-Singer es una marca comercial registrada de Hansa Medical Products, Inc., InHealth es una marca comercial registrada de International Healthcare Technologies, División de Helix Medical, Inc., / Blom-Singer est une marque déposée de Hansa Medical Products, Inc., InHealth est une marque déposée de International Healthcare Technologies, filiale de Helix Medical, Inc., / Blom-Singer è un marchio depositato di Hansa Medical Products, Inc., InHealth è un marchio depositato di International Healthcare Technologies, una divisione di Helix Medical, Inc., / Blom-Singer is een gedeponeerd handelsmerk van Hansa Medical Products, Inc., InHealth is een gedeponeerd handelsmerk van International Healthcare Technologies, een divisie van Helix Medical, Inc., / Blom-Singer é uma marca registada da Hansa Medical Products, Inc., InHealth é uma marca registada da International Healthcare Technologies, uma Divisão da Helix Medical, Inc.