

# Styles and Features



Presented by Eric D. Blom, Ph.D.

The information in this chart has been compiled by Dr. Blom based upon his extensive clinical experience. Refer to product data sheet for complete information.

**Blom-Singer®**

voice restoration systems

*The Voice of Experience  
Since 1978*



## Voice Prostheses

### FEATURES

### STYLES

**ADVANTAGE™**  
Indwelling Voice Prosthesis  
DIAMETERS AVAILABLE  
16 Fr. & 20 Fr.

**CLASSIC**  
Indwelling Voice Prosthesis  
DIAMETERS AVAILABLE  
16 Fr. & 20 Fr.

**Low Pressure Voice Prosthesis**  
DIAMETERS AVAILABLE  
16 Fr. & 20 Fr.

**Duckbill Voice Prosthesis**  
DIAMETERS AVAILABLE  
16 Fr.

Patient can remove/insert voice prosthesis.

Patient cannot remove/insert voice prosthesis.

Employs atraumatic Blom-Singer Gel Cap Insertion System.

Small diameter prosthesis (16 Fr.) for use in questionably healthy TE wall tissue.

Anti-fungal additive (silver oxide) incorporated in valve as a material preservative.

Voice prosthesis usually requires no neck strap or adhesive tape.


Esophageal flange provides maximal retention.

Valve (20 Fr.) offers minimal resistance to speech airflow.

Voice prosthesis (20 Fr.) is available in one millimeter (mm) incremental lengths.

Optional increased valve resistance to prevent air ingestion to stomach.

Optional large esophageal flange (20 Fr.) to prevent leakage around prosthesis.



(20 Fr.)